

<b>CAMERA</b>	<b>HEAD SIZE</b>	<b>RECOMMENDED NOMINAL ACTIVITY</b>	<b><u>C-THRU</u> Active Area</b>	<b><u>FEATHERLITE</u> Active Area</b>	<b><u>STANDARD</u> Active Area</b>	<b><u>PERFLEXION</u> Active Area</b>	<b>CAMERA SPECIFIC</b>
<b>Digirad</b> 2020tc	8" x 8"	10 mCi	C-THRU <b>CTRF10012</b> 8.9" x 8.9"	FEATHERLITE <b>MED3722</b> 11" x 11"	N/A	PERFLEXION <b>PF09R</b> 9" x 9"	N/A
<b>Digirad</b> Cardius Series	8" x 6"	10 mCi	C-THRU <b>CTRF10012</b> 8.9" x 8.9"	FEATHERLITE <b>MED3722</b> 11" x 11"	N/A	N/A	N/A
<b>Digirad</b> SPECT Pak	9" x 9"	10 mCi	C-THRU <b>CTRF10012</b> 8.9" x 8.9"	FEATHERLITE <b>MED3722</b> 11" x 11"	N/A	PERFLEXION <b>PF09R</b> 9" x 9"	N/A
<b>Dilon</b> 6800	8" x 6"	5 mCi	C-THRU <b>CTRF10012</b> 8.9" x 8.9"	FEATHERLITE <b>MED3715</b> 8.22" x 14.25"	N/A	N/A	N/A
<b>Elscint</b> Cardial	15.75" x 10"	10 mCi Max	C-THRU <b>CTRF10017</b> 16" x 10"	FEATHERLITE <b>MED3713</b> 18" x 14"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A
<b>Elscint</b> Helix	21.25" x 15.75"	10 mCi Max	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Elscint</b> SP4	15.75"	10 mCi Max	C-THRU <b>CTRF10050</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES392</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
<b>Elscint</b> SP6	21.25" x 15.75"	10 mCi Max	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Elscint</b> Varicam	21.25" x 15.75"	10 mCi Max	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>GE</b> 300 AM	10"	3 mCi	N/A	FEATHERLITE <b>MED3700</b> 14"	STANDARD <b>NES298</b> 14"	N/A	N/A
<b>GE</b> 5000 Series	24"	10 mCi	N/A	FEATHERLITE <b>MED3708</b> 24"	N/A	N/A	N/A

<b>GE</b> Discovery 530	14" x 7.5"	10 mCi	C-THRU <b>CTRF10017</b> 16" x 10"	FEATHERLITE <b>MED3727</b> 16" x 10"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A
<b>GE</b> Discovery VH	20" x 14"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>GE</b> DSTi/Dsi	14" x 13"	10 mCi	N/A	FEATHERLITE <b>MED3713</b> 18" x 14"	STANDARD <b>NES8480</b> 18" x 14"	N/A	N/A
<b>GE</b> Hawkeye	20" x 14"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>GE</b> Infinia	21.25" x 15.75"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>GE</b> Maxi 2	15"	5 mCi	C-THRU <b>CTRF10051</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES391</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
<b>GE</b> Maxi 37	15"	5 mCi	C-THRU <b>CTRF10050</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES391</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
<b>GE</b> Maxxus	21" x 16"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>GE</b> Millenium MPS	14" x 14"	10 mCi	N/A	FEATHERLITE <b>MED3725</b> 14" x 14"	STANDARD <b>NES8480</b> 18" x 14"	N/A	N/A
<b>GE</b> Millenium MPR/VG	21.75" x 15.75"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>GE</b> Millenium MG	20" x 14"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>GE</b> Myosight	20" x 14"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A

<b>GE</b> Neurocam	7.9" x 6.7"	10 mCi	C-THRU <b>CTRF10012</b> 8.9" x 8.9"	FEATHERLITE <b>MED3722</b> 11" x 11"	N/A	PERFLEXION <b>PF09R</b> 9" x 9"	N/A
<b>GE</b> Optima	14" x 9"	7.5 mCi	C-THRU <b>CTRF10017</b> 16" x 10"	FEATHERLITE <b>MED3727</b> 16" x 10"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A
<b>GE</b> Starcam/ACT	15.4"	10 mCi	C-THRU <b>CTRF10050</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES392</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
<b>GE</b> Starcam/XCT	15.4"	10 mCi	C-THRU <b>CTRF10050</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES392</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
<b>GE</b> Starcam/XRT	21" x 16"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>GE</b> Ventri	14.6" x 7.4"	10 mCi	C-THRU <b>CTRF10017</b> 16" x 10"	FEATHERLITE <b>MED3727</b> 16" x 10"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A
<b>Hitachi</b> 1024C	13.5"	5 mCi	N/A	FEATHERLITE <b>MED3700</b> 14"	STANDARD <b>NES298</b> 14"	N/A	N/A
<b>Hitachi</b> 1024R/1024RDT	19.7" x 14.2"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Hitachi</b> 150, 250DSP	20" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Hitachi</b> 260DSP	22" x 16"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>Hitachi</b> NeuroSpect	8.7" x 6.7"	10 mCi	C-THRU <b>CTRF10012</b> 8.9" x 8.9"	FEATHERLITE <b>MED3722</b> 11" x 11"	N/A	PERFLEXION <b>PF09R</b> 9" x 9"	N/A
<b>IS2</b> DCC	15" x 10.2"	10 mCi	C-THRU <b>CTRF10017</b> 16" x 10"	FEATHERLITE <b>MED3727</b> 16" x 10"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A

<b>IS2</b> Pulse CDC	15" x 10.2"	10 mCi	C-THRU <b>CTRF10017</b> 16" x 10"	FEATHERLITE <b>MED3727</b> 16" x 10"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A
<b>IS2</b> SR	20" x 14.8"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>IS2</b> SC	15"	10 mCi	C-THRU <b>CTRF10050</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES392</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
Mediso CardioSpect D90	14.5" x 9"	10 mCi	C-THRU <b>CTRF13770A</b> 15.2" x 9.1"	FEATHERLITE <b>MED3727</b> 16" x 10"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A
Mediso CardioSpect SC	20.8" x 15.3"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
Mediso CardioSpect SR	20.8" x 15.3"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
Mediso CardioSpect VMAX	20.8" x 15.3"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
Mediso NeuroSpect Quad	9" x 8"	10 mCi	C-THRU <b>CTRF13770A</b> 15.2" x 9.1"	FEATHERLITE <b>MED3722</b> 11" x 11"	N/A	PERFLEXION <b>PF09R</b> 9" x 9"	N/A
<b>Park Medical</b> Isocam I & II	22.3" x 16.5"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (ADAC)</b> ARC 3000	15"	10 mCi	C-THRU <b>CTRF10050</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES392</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
<b>Philips (ADAC)</b> Argus	20" x 15"	10 mCi Max	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (ADAC)</b> Brightview	20.25" x 16"	10 mCi Max	C-THRU <b>CTRF10021</b> 24.4" x 17.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8610</b> 24" x 17.5"	PERFLEXION <b>PF24RW</b> 23.9" x 16.4"	N/A

<b>Philips (ADAC)</b> Cardial	20" x 15"	10 mCi Max	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (ADAC)</b> Cardio 60	20" x 15"	10 mCi Max	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (ADAC)</b> Cardio MD	9.2" x 15.4"	10 mCi Max	C-THRU <b>CTRF10017</b> 16" x 10"	FEATHERLITE <b>MED3727</b> 16" x 10"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A
<b>Philips (ADAC)</b> Cirrus	15"	10 mCi	C-THRU <b>CTRF10050</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES392</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
<b>Philips (ADAC)</b> Forte	20" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (ADAC)</b> Genesys	20" x 15"	10 mCi Max	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (ADAC)</b> Meridian	20" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (ADAC)</b> Precedence	20" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (ADAC)</b> Skylight	20" x 14"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (ADAC)</b> Solus	20" x 15"	10 mCi Max	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (ADAC)</b> Vertex	20" x 15"	10 mCi Max	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (Picker)</b> 411	11"	5 mCi	N/A	FEATHERLITE <b>MED3700</b> 14"	STANDARD <b>NES297</b> 14"	N/A	N/A

<b>Philips (Picker)</b> 412	12"	5 mCi	N/A	FEATHERLITE <b>MED3700</b> 14"	STANDARD <b>NES297</b> 14"	N/A	N/A
<b>Philips (Picker)</b> 415	15"	5 mCi	C-THRU <b>CTRF10051</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES391</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
<b>Philips (Picker)</b> Axis/Irix	21" x 15.5"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>Philips (Picker)</b> Dynamo	10"	3 mCi	N/A	FEATHERLITE <b>MED3700</b> 14"	STANDARD <b>NES297</b> 14"	N/A	N/A
<b>Philips (Picker)</b> PRISM 1000	20" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (Picker)</b> PRISM 2000	20" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Philips (Picker)</b> PRISM 3000	15.7" x 9.4"	10 mCi	C-THRU <b>CTRF10017</b> 16" x 10"	FEATHERLITE <b>MED3727</b> 16" x 10"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A
<b>Philips (Picker)</b> SX300	14" x 14"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3725</b> 14" x 14"	STANDARD <b>NES8480</b> 18" x 14"	N/A	N/A
<b>Siemens</b> Bodyscan	23.5" x 15.75"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Siemens</b> c.cam	14" x 8.4"	10 mCi	C-THRU <b>CTRF10017</b> 16" x 10"	FEATHERLITE <b>MED3727</b> 16" x 10"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A
<b>Siemens</b> Diacam	21.25" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>Siemens</b> e.cam	21.25" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"

<b>Siemens</b> MultiSpect 2	21.25" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>Siemens</b> MultiSpect 3	16" x 12"	10 mCi	N/A	FEATHERLITE <b>MED3713</b> 18" x 14"	STANDARD <b>NES8480</b> 18" x 14"	N/A	N/A
<b>Siemens</b> Orbiter 3700, 7500	15.25"	10 mCi	C-THRU <b>CTRF10050</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES392</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
<b>Siemens</b> Orbiter 3700, 7500	21.25" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>Siemens</b> Symbia	21.23" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>SMV (Sopha)</b> BodyTrac	21.3" x 15.75"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>SMV (Sopha)</b> DS7	15.75"	10 mCi	C-THRU <b>CTRF10050</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES392</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
<b>SMV (Sopha)</b> DST	15.75" x 11.8"	10 mCi	N/A	FEATHERLITE <b>MED3713</b> 18" x 14"	STANDARD <b>NES8480</b> 18" x 14"	N/A	N/A
<b>SMV (Sopha)</b> DSX	21.3" x 15.75"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>SMV (Sopha)</b> DSTXL	21.3" x 15.75"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>SMV (Sopha)</b> Vision FX	20" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>Technicare</b> Omega 500	20" x 14.5"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A

<b>Technicare</b> Gemini 600	15"	10 mCi	C-THRU <b>CTRF10050</b> 18.5"	FEATHERLITE <b>MED3706</b> 18.5"	STANDARD <b>NES392</b> 18.5"	PERFLEXION <b>PF18C</b> 18.5"	N/A
<b>Technicare</b> Gemini 700	20" x 14.5"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Toshiba</b> GCA 601A	13.75"	5 mCi	N/A	FEATHERLITE <b>MED3700</b> 14"	STANDARD <b>NES297</b> 14"	N/A	N/A
<b>Toshiba</b> GCA 602/602A	13.75"	5 mCi	N/A	FEATHERLITE <b>MED3700</b> 14"	STANDARD <b>NES297</b> 14"	N/A	N/A
<b>Toshiba</b> 7100A/7200	21.5" x 19"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>Toshiba</b> GCA 901/901A	20.1" x 14.5"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Toshiba</b> 901WB	20.1" x 14.6"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Toshiba</b> GCA9300	16.5" x 8.3"	10 mCi	C-THRU <b>CTRF10017</b> 16" x 10"	FEATHERLITE <b>MED3727</b> 16" x 10"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A
<b>Toshiba</b> t.cam	21.25" x 15"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	SPECIFIC <b>JRT-PZ5</b> 22.2" x 16.7"
<b>Trionix</b> Biad 20/XLT20	20" x 14.5"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Trionix</b> Biad 24/XLT24	24" x 14.5"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A
<b>Trionix</b> Monad	24" x 14.5"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A

<b>Trionix</b> Triad XLT	15.75" x 8.7"	10 mCi	C-THRU <b>CTRF10017</b> 16" x 10"	FEATHERLITE <b>MED3727</b> 16" x 10"	STANDARD <b>NES8510</b> 16" x 10.5"	PERFLEXION <b>PF16R</b> 16" x 10.5"	N/A
<b>Trionix</b> Triad XLT20	20" x 14.5"	10 mCi	C-THRU <b>CTRF10000</b> 24.4" x 16.5"	FEATHERLITE <b>MED3709</b> 24" x 16.5"	STANDARD <b>NES8400</b> 23.9" x 16.4"	PERFLEXION <b>PF24R</b> 23.9" x 16.4"	N/A

|

|

|



|

|

|



|