

XENON CONVENIENCE KITS™

Featuring an air-cushioned face mask, clear ultra-flex tubing which extends from 6" to 24" and mouthpiece kits with direct dose administration adapter for leak-proof xenon delivery


Biodex Xenon Convenience Kits bring together all the components needed to complete a single xenon study. Xenon can be administered via direct dose administration, that conveniently luer locks without the use of a needle; or by injection port, which requires a syringe needle.

Easy to use, simple to order and disposable, Biodex Xenon Convenience Kits are more than a bargain, they are a sensible

and time-saving investment. In fact, our complete kits can prolong the life of your xenon charcoal traps by ensuring that the pre-filled Drierite cartridges are fresh for every scan while the pre-filled Soda-Lime (CO₂ absorber) cartridges eliminate possible breakdown of granules that can lead to a clogged system and hinder patient air flow.



FIND OUT HOW
TO QUALIFY FOR A
50% DISCOUNT*
ON A
CHARCOAL TRAP



XENON CONVENIENCE KITS™

AIR-CUSHIONED™ Face Mask:

- 132-680** Convenience Kit, Face Mask ☉
Includes: Bacteria filter
- 132-681** Convenience Kit, Face Mask ☉
Includes: Bacteria filter and ultra-flex tubing (shown)
- 132-781** Convenience Kit, Face Mask ☉
Includes: Bacteria filter, ultra-flex tubing, Drierite & Soda-Lime cartridges

XENON CONVENIENCE KITS™

AIR-CUSHIONED™ Face Mask w/Injection Port:

- 132-690** Convenience Kit, Face Mask with Injection Port ☉
Includes: Bacteria filter (shown)
- 132-691** Convenience Kit, Face Mask with Injection Port ☉
Includes: Bacteria filter and ultra-flex tubing
- 132-784** Convenience Kit, Face Mask with Injection Port ☉
Includes: Bacteria filter, ultra-flex tubing, Drierite & Soda-Lime cartridges

CALL FOR FREE SAMPLE*

☉ DENOTES SINGLE USE ONLY

* Offer applies to the United States only.

XENON CONVENIENCE KITS™

Direct Dose Administration for leak-proof xenon delivery



SECTION 7

XENON CONVENIENCE KITS™

AIRCUSHIONED™ Face Mask
w/Direct Dose Administration:

- 132-699** Convenience Kit, Face Mask with Direct Dose ©
Includes: Bacteria filter
- 132-692** Convenience Kit, Face Mask with Direct Dose ©
Includes: Bacteria filter and ultra-flex tubing
- 132-793** Convenience Kit, Face Mask with Direct Dose ©
Includes: Bacteria filter, ultra-flex tubing, Drierite & Soda-Lime cartridges (shown)



XENON CONVENIENCE KITS™

with Mouthpiece and Direct Dose Administration Adapter:

- 132-770** Convenience Kit, Mouthpiece ©
Includes: Bacteria filter and nose clip
- 132-774** Convenience Kit, Mouthpiece and Administration Adapter ©
Includes: Bacteria filter, nose clip and ultra-flex tubing (shown)
- 132-771** Convenience Kit, Mouthpiece and Administration Adapter ©
Includes: Bacteria filter, nose clip, ultra-flex tubing, Drierite & Soda-Lime cartridges

XENON CONVENIENCE KIT™

AIRCUSHIONED™ Face Mask and Direct Dose Administration Adapter:

- 134-772** Convenience Kit, Face Mask and Administration Adapter ©
Includes: Bacteria filter and ultra-flex tubing

CALL FOR FREE SAMPLE*

© DENOTES SINGLE USE ONLY


* Offer applies to the United States only.

This product is available through:

JRT Associates

800-221-0111

AIR-CUSHIONED FACE MASK



The air-cushioned face mask has a pre-filled air-cushion (medium inflation pressure) that molds to the contour of the patient's face providing a tight, leak-proof and comfortable seal. Cushion pressure is easily adjusted by inserting a standard syringe without needle in the two-way valve located on the mask bottom. The flexibility of cushion pressure allows optimum surface contact for every patient. Transparency of the entire mask allows continuous visual identification of patient's vital signs.

AIR-CUSHIONED™ Face Mask without Injection Port:
132-685 Face Mask without Injection Port, Adult ☉



AIR-CUSHIONED™ Face Mask with Injection Port:
132-695 Face Mask with Injection Port, Adult ☉

AIR-CUSHIONED FACE MASK

with Direct Dose Administration




With direct dose administration

In addition to the patient comfort provided by the soft, pliable air cushion of the disposable face mask, the direct dose feature reduces patient anxiety because no needle is used. The syringe luer locks tightly to the direct dose tube, away from the patient's face. The one-way valve prevents the gas from reversing through the mask luer connector, allowing the technologist to disconnect immediately after injecting.

AIR-CUSHIONED™ Face Mask with Direct Dose Administration:
132-698 Face Mask with Direct Dose Administration, Adult ☉

FACE MASK HARNESS



This traditional Face Mask Harness is made of soft rubber and can be adjusted for any head size. The square rubber base is comfortably positioned on the back of the patient's head and the harness tails are brought around to the front. The small holes fit snugly on the hook ring and hold the mask firmly in place.

139-677 Face Mask Harness, Adult

SECTION 7

☉ DENOTES SINGLE USE ONLY

DISPOSABLE MOUTHPIECE



This mouthpiece has a soft flange, formed tooth-grip protrusions, and fits snugly between the patient's lips and teeth to form a tight seal.

130-550 Mouthpiece, Disposable ©

TRU-FIT DISPOSABLE MOUTHPIECE



Contoured to the natural shape of the mouth, this mouthpiece can be held gently, but securely — eliminating the need to clamp down or the possibility of “popping” out of the patient's mouth and interrupting the study. The soft plastic material is transparent so that any obstruction can be immediately detected by the technologist. Pleasantly scented

with vanilla, this mouthpiece makes the procedure a little more comfortable for the patient.

130-551 Mouthpiece, Disposable, Vanilla Scented ©


DISPOSABLE NOSE CLIP



Economical, this nose clamp is all plastic for single patient use. Used with any mouthpiece.

130-100 Nose Clip, Disposable, 100/pkg ©

DISPOSABLE BACTERIA FILTER



Bacteria filters are used to reduce the possibility of cross contamination. The single-use filter is placed in line between a delivery tube and disposable mouthpiece or face mask. Electrostatically charged filter media is 99.9% effective in bacteria/virus retention while maintaining low breathing resistance. Two filters can be piggy-backed together for potential high-risk studies.

132-750 Bacteria Filter ©

BLANKET ORDERS

You may qualify to receive a 50% discount on a replacement Charcoal Trap for your Pulmonex just by placing a Blanket Order.*

- There is no binding contract. You can cancel at any time without penalty or change the quantities throughout the life of the blanket.
- You set the quantities and delivery intervals — based on your department needs.
- A single purchase order establishes a tailored blanket for up to one year.
- A Blanket Order exempts you from price increases on blanketed items.
- Disposables are there when you need them. No more rush orders, expensive overnight shipping charges or equipment downtime.
- Maximize your storage space — by using ours.

**Offer applies to United States only.*

© DENOTES SINGLE USE ONLY

* Offer applies to the United States only.

This product is available through:

JRT Associates

800-221-0111

CORRUGATED TUBING



Corrugated Tubing is scored and capped at 6" intervals for easy cutting and firm attachment to fittings.

139-680 Tubing, Corrugated, 100 ft/roll ☉

TUBING SPLICE



Tubing Splice is used to connect sections of 22 mm corrugated tubing to individual mask and filter.

130-639 Tubing Splice, 22 mm Male, 10/pkg ☉

"Y" CONNECTOR



The "Y" Connector has a one-way valve and plug. It is used to connect two 22 mm tubes to a mouthpiece, face mask or bacteria filter.

185-302 "Y" Connector, luer plug ☉

139-102 "Y" Connector, solid plug ☉

☉ DENOTES SINGLE USE ONLY

DRIERITE

Moisture absorber - new blue and white blend



Drierite serves as a moisture trap for the air going into the charcoal trap of the Pulmonex Xenon System. Drierite is blue and white when dry and turns pink when it has absorbed maximum moisture. Available in two convenient sizes, Drierite is packaged in an air-tight container with a wide-access mouth.

139-101 Drierite, 1 lb

139-104 Drierite, 4.5 lb

132-772 Drierite, Disposable Cartridge ☉

SODA-LIME

CO₂ absorber



Soda-Lime, a highly efficient CO₂ absorber, prevents patients from re-breathing carbon dioxide and subsequent acidosis. Available in two convenient sizes, Soda-Lime is packaged in an air-tight container with a wide-access mouth.

130-019 Soda-Lime, 1 lb

130-020 Soda-Lime, 4.5 lb

132-773 Soda-Lime, Disposable Cartridge ☉